


REDUCE EPO USAGE with Nephros Dual Stage Ultrafilter



Achieve Significant Benefits Beyond Using
Standard Dialysate Ultrafilters Alone

OBSERVED DECREASE OF EPO USAGE: 25 - 50%

**Erythropoietin (EPO) Use Reduction
with Nephros Dual Stage Ultrafilter (DSU)**



Clinic 1

- Baseline with DiaClear™
- 6 Month Study
- 52% EPO Use Reduction
- $p = 0.001$

Clinic 2

- Baseline without ultrafilter
- 3 Month Study
- 26% EPO Use Reduction
- $p = 0.0002$

Clinic 3

- Baseline with DiaSafe™
- 3 Month Study
- 35% EPO Use Reduction
- $p = 0.003$

* DSU ultrafilters installed on incoming water line to dialysis machine