

**DSU FILTER**

**EAST COAST DIALYSIS UNIT - JULY 2010 - JANUARY 2011**

Pt #	EPO DOSAGE(1000 UNITS)						2011	HGB							2011
	JULY	AUG	SEPT	OCT	NOV	DEC	JAN	JULY	AUG	SEPT	OCT	NOV	DEC	JAN	
1-085	22000	22000	16500	8800	7700	7700	7000	11.3	11.7	12.6	12.9	11.9	11.3	12.9	
7-090	8800	11000	11000	11000	11000	8000	8000	12.5	10.7	12.1	10.9	13.4	12.8	12.6	
8-111	11000	7700	6600	5500	5500	4000	4000	12.7	12.0	11.9	12.8	13.2	12.0	11.2	
9-077	33000	17600	22000	20900	20600	17600	17600	12.3	11.7	10.1	11.5	11.0	11.3	9.0	
10-043	12100	13200	12100	9900	5500	7000	7000	11.3	10.7	11.6	11.9	11.9	10.8	8.8	
11-023	11000	15400	18700	15400	15400	7000	7000	10.6	10.3	10.2	11.7	11.7	11.5	12.3	
13-125	16500	16500	16500	12100	12100	15400	15400	12.1	10.5	11.1	12.2	9.4	10.2	9.3	
14-082	11000	19800	17600	13200	13200	4000	4000	10.0	9.8	13.3	11.1	11.2	10.8	10.7	
15-084	8800	4400	2200	2200	2200	2000	2000	12.8	11.9	12.4	11.6	11.5	10.3	10.5	
16-024	6600	6600	4400	2200	2200	3000	3000	10.0	11.1	11.3	12.3	11.6	10.8	10.0	
18-040	7700	4400	2200	2200	2200	2000	2000	12.7	12.8	12.5	12.4	11.6	10.4	10.1	
20-101	8800	5500	4400	2200	2200	2200	6000	12.0	10.7	12.8	12.5	12.2	9.9	11.9	
21-103	15400	11000	6600	0	6600	6000	6000	11.9	11.2	14.1	12.0	9.1	9.7	10.4	
31-049	11000	15400	12100	12000	8800	7700	5500	10.0	10.2	12.5	13.1	12.0	12.3	11.8	
34-004	5500	5500	5500	3300	3300	5000	5000	10.7	11.2	11.5	12.3	10.8	8.7	9.7	
37-054	5500	5500	5500	4400	4400	3000	3000	10.7	11.4	13.0	11.3	11.7	11.7	12.8	
39-110	7700	7700	5500	5500	4400	4400	2000	12.3	12.7	11.5	13.3	13.1	13.2	11.3	
40-005	15400	11000	6600	8800	8800	6000	6000	12.4	13.3	10.6	10.4	10.5	11.5	11.2	
45-098	4400	6600	5500	5500	5500	5500	5500	11.1	11.4	10.9	11.7	13.5	12.5	12.4	
48-057	4400	4400	4400	4400	4400			11.8	11.3	11.9	11.1	8.7			
49-114	5500	4400	3300	2200	2200	2200	2000	10.6	12.2	12.2	11.8	11.7	12.6	12.4	
<b>TOTALS</b>	<b>232100</b>	<b>215600</b>	<b>189200</b>	<b>151700</b>	<b>148200</b>	<b>119700</b>	<b>118000</b>	241.8	238.8	250.1	250.8	241.7	224.3	221.3	
<b>AVG</b>	<b>11052</b>	<b>10267</b>	<b>9010</b>	<b>7224</b>	<b>7057</b>	<b>5985</b>	<b>5900</b>	<b>11.5</b>	<b>11.4</b>	<b>11.9</b>	<b>11.9</b>	<b>11.5</b>	<b>11.2</b>	<b>11.1</b>	
							<b>47%</b>							<b>4%</b>	
							<b>DECREASE</b>							<b>REDUCTION</b>	

Note: Fresenius Diasafe filters were used prior to installation of the Nephros DSU in July 2010.